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Kit of reagents for the execution of Voges-Proskauer test

DESCRIPTION

V.P. TEST EP Reagent is a kit containing the reagents ALPHA NAPHTHOL and KOH 40% necessary for the execution of the Voges-Proskauer test, used for identification of microorganisms performing 2-3 butilneglicole fermentation.

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Each package contains 1 INSTRUCTION SHEET, 2 DROPPERS and the follow reagents:

ALPHA NAPHTHOL	KOH 40%
----------------	---------

METHOD PRINCIPLE

Voges-Proskauer reaction distinguishes organisms that ferment glucose to acids from those that ferment to acetoin (2-3 butilenglicole fermentation). After 48 hours incubation ALPHA NAPHTHOL and KOH 40% reagents are added to sample. Positive reaction to Voges-Proskauer test is highlighted by the formation of a red coloration obtained from acetyl-methyl-carbinol oxidation in the presence of alkali and oxygen.

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Reagents	Composition	Contents
ALPHA NAPHTHOL	Alpha Naphthol 6 g, Ethyl Alcohol 100 ml	5 Vials X 10 ml
KOH 40%	Potassium Hydroxide 40 g, Distilled Water 100 ml	5 Vials X 10 ml

TEST PROCEDURE AND RESULTS INTERPRETATION

The VP test is performed on cultures of 48 hours of incubation prepared dissolving a pure culture in a tube of MR-VP BROTH (ref. 610032 or 620032) properly prepared, or in a tube of MR-VP MEDIUM (ref. 20149). Add 3 drops of ALPHA NAPHTHOL and 2 drops of KOH 40% solution to 3 ml of broth culture. The positive test is given by the development of a bright red colour 15 minutes after addition of reagents.

QUALITY CONTROL FOR THE USER

Visual inspection

KOH 40% is a colourless clear solution;

ALPHA NAPHTHOL is a beige colour clear solution.

Microbiology control

Each batch of V.P. TEST EP Reagent is subjected to quality control using a bacterial culture of *Escherichia coli* ATCC 25922 as negative control and one of *Enterobacter cloacae* ATCC 13047 as positive control.

PRECAUTIONS

The product V.P. TEST EP Reagent, is classified as hazardous under current legislation; for its use it is recommend to consult the MSDS. V.P. TEST EP Reagent is a disposable device to be used only for in vitro diagnostic use, is intended for professional use and should be used in the laboratory by properly trained operator, with approved methods of aseptic and safety against pathogens.

STORAGE AND TRANSPORT CONDITIONS

Store V.P. TEST EP Reagent at 2-8°C away from light in its original package, until the expiry date shown on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATION OF THE USED MATERIAL

After use, V.P. TEST EP Reagent and the material that comes in contact with the sample must be decontaminated and disposed in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infectious materials.

BIBLIOGRAPHY

1. Edwin H. Lenette: Manual of Clinical Microbiology (1995).
2. Blazevic, D.J., and Ederer, G.M.: Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons, 1975.

PRESENTATION

Product	REF	
V.P. TEST EP Reagent	80281	100 tests

TABLE OF SYMBOLS

LOT Batch code	IVD In vitro Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse